Automated processing of solid tissues into single cells or nuclei for genomics and cell biology applications with the Singulator[™] 100 and 200 systems.



S. Jovanovich*, J. Bashkin, K. Chear, N. Khan, K. Kumbagowdana, B. Leisz, R. McIntosh, D. Meyer, N. Pereira S2 Genomics, Livermore, CA.

Single-cell sequencing is revealing the next level of complexity in biological systems. Single-cell sequencing workflows for multiple applications require reproducible generation of high-quality single-cell or nuclei suspensions.

S2 Genomics developed and Singulator Cartridge commercialized the patented Singulator[™] 100 System to automate the dissociation of solid tissues into single cell or nuclei suspensions in single-use cartridges using coupled enzymatic 100 Singulator™ Figure System, or chemical dissociation and mechanical Singulator including disruption. S2 Genomics is now developing instrument, the the Singulator 200 System[™], with higher single-use cartridge, chiller for nuclei sample throughput, improved temperature reagents, and Single-Shot Mechanism for and force control, and other advancements. cell preparation reagents. Validation data will be presented for processing <1 to 400 mg samples of fresh 180,000 160,000 mouse and human solid tissues into single-cell 140,000 suspensions with high viabilities and yields. In ີ ພິ 120,000 **a** 100,000 addition, the Singulator can process fresh or 80,000 frozen mouse and human tissue samples into 60,000 high-quality nuclei in ~5 min. 40,000 20,000 Singulator-prepared samples have been analyzed by single-cell RNA sequencing Spleen Kidney Brain Lung Gut Liver (1:10) (scRNA-seq), single-nuclei RNA sequencing **Tissue Type** (snRNA-seq), and assay for transposase-100% accessible chromatin using sequencing (ATAC-90% 80% seq). 70% 60% Systems can Singulator The 50% 40% automatically process fresh tissue 30% 20% samples into single cell suspensions, 10% while nuclei can be isolated from Kidney Spleer Brain Liver (1:10) **Tissue Type** fresh, frozen, or OCT preserved tissue. Figure 1 shows the Singulator Figure 3. Automated dissociation of 100 for processing a single sample. mouse tissues into single cells. Tissue-Figure 2 shows the Singulator 200 high viabilities (bottom). capable of processing two samples 350,000 with random access. The Singulator 300,000 has SingleSense™ real-time 200 250,000 pressure and temperature control for ß 200,000 precise, gentle dissociation. **..** 150,000 The Singulator platforms are open 100,000 with respect to reagents: use S2 50,000 reagents preferred your or formulations from any source. (1:10)Enhanced software let's you create **Tissue Type** dissociation of Figure 4. Automated and share the right process for your mouse tissues into nuclei. Tissue-specific tissue. yields are obtained in ~five minutes.





Figure from data shows processing fresh mouse samples into cell suspensions using the Singulator 100 System with standard cartridges. Samples that are difficult to process into single cells, such as frozen samples, or human prostate, heart, or muscle, can be quickly and routinely processed into nuclei suspensions for ATAC-Seq or snRNA-Seq analysis (Figure 4).

Figure 2. Singulator[™] 200 System for single-use processing two samples with chiller for cartridges, integrated Single-Shot nuclei reagents, and dual



delivering Mechanisms for preparation reagents.



Figure 5. Dissociation of small samples with the NIC+ cartridges. Yield of mouse kidney nuclei per mg vs. sample size is shown for standard (yellow) and NIC+ (green) cartridges.

recently introduced NIC+ The cartridges produce nuclei efficiently with samples as small as 1 mg, including biopsy samples. Figure 5 illustrates the improved efficiency of and yields of the NIC+ cartridges for kidney nuclei isolations. mouse Similar results are obtained from frozen samples and other tissues.

The Singulator Systems are being applied for a wide range of studies, from basic biology to translational research. Figure 6 shows scRNA-Seq from infiltrating tumor results leukocytes (TILs) purified from a human lung tumor with magnetic beads.





cell

Figure 6. scRNA-Seq of TILS. Human lung samples were dissociated with the Singulator 100 with the cell isolation cartridges, and TILs purified by magnetic separation (StemCells) before scRNA-Seq library construction with a Chromium. 7,676 TIL cells were sequenced with analysis by Azimuth using the PBMC data set.

Summary: The Singulator 100 System automates processing fresh or frozen tissue into filtered suspensions of single cells or nuclei. Samples down to 1 mg can now be routinely processed for nuclei analysis.

*stevanj@s2genomics.com

Supported in part by National Institutes of Health award number R41 HG010129-01 to Dr. Jovanovich.