

USER GUIDE

Curio Seeker 10x10 Spatial Mapping Kit

For fresh frozen tissues

Table of Contents

Version summary	3
Overview	4
Estimated workflow timing	6
Materials	7
Curio Seeker 10x10 Kit components (4 reactions)	7
Variations of the Curio Seeker Dual Indexing Primer Kit: v1 and v2	8
Curio Seeker Dual Indexing Primer Kit v1 components	8
Curio Seeker Dual Indexing Primer Kit v2 components	9
Curio Seeker reaction chamber pouch components (4 reactions)	9
Curio Seeker 10x10 adapter	10
Required reagents not included	10
Equipment	11
Consumables	12
Tissue requirements and recommendations for assessing tissue quality	12
Tips and techniques	13
Protocol	14
A. Preparation before starting	14
B. Tissue sectioning and hybridization to Curio Seeker tile	14
C. Reverse transcription	17
D. Tissue clearing and Curio Seeker bead resuspension	18
E. Second strand synthesis	21
F. cDNA amplification	22
G. Purification and quantification	24
H. Tagmentation (with Nextera XT Library Prep Kit)	26
I. Library cleanup and QC	27
J. Sequencing	30
K. Informatics	32
Troubleshooting	33
Appendix A: Indexing PCR with Curio Seeker Dual Indexing Primer Kit (K003)	36

Version summary

Revision date	Version	Major changes
Oct 2023	1.1	<ul style="list-style-type: none">• Edited Appendix with minor changes• Add PCR tube magnetic rack• Correction made to warning on page 27
Feb 2024	1.2	<ul style="list-style-type: none">• Added static electricity warning• Added storage temperature recommendations for master mixes• Corrected Curio Seeker Dual Indexing Primer Kit v1 instructions• RIN score addendum• Added CryoCube guidelines in section B• Increased hybridization time from 15 minutes to 30 minutes• Added buffer temperature guidelines for the bead dissociation step in section D
May 2024	1.2b	<ul style="list-style-type: none">• Minor edition of schematics

Overview

This user manual describes the 10x10 Curio Seeker Spatial Mapping Kit protocol, hereby referred to as the Curio Seeker workflow to generate high-quality, Illumina sequencing-ready libraries from fresh frozen tissues to obtain high resolution spatial transcriptomic information of a sample. Once tissues are sectioned and placed onto the Curio Seeker tile (hereby referred to as 'tile') (Figure 1), a glass substrate containing a monolayer of uniquely DNA-barcoded microparticles (referred to as 'beads' in this document), the rest of the workflow can be completed in under eight hours (Figure 2), with multiple safe stopping points. The Curio Seeker workflow (Figure 3) starts with hybridization of RNA to the beads on the tile, followed by reverse transcription. A tissue clearing step is performed to digest the tissue and release the beads from the glass into solution. Next, second strand synthesis is performed by semi-random priming followed by cDNA amplification. Finally, the Nextera XT DNA Sample Preparation Kit is used to generate Illumina sequencing platform-compatible libraries.




Figure 1. Curio Seeker tile

Curio Seeker workflow




Figure 3. Curio Seeker workflow. Please refer to supplementary sections at the end of this protocol for more detailed sequence information.

Estimated workflow timing



















Figure 3. Estimated workflow timing

Materials











Curio Seeker 10x10 Kit components (4 reactions)

Component	Part number
Pouch 1 (4°C)	LTB001
10x10 Seeker Indexed Tile	SL001
Box 1 (RT)	K001
 Hyb Buffer	B005
 TC Buffer	B004
 Nuclease-free Water	B001
 Bead Wash Buffer	B003
Box 2 (-20°C)	K002
 RNase Inhibitor	E001
 RT/SS Buffer	B006
 dNTP	N001
 RT Enzyme	E003
 TC Enzyme	E002
 SS Primer	P002
 SS Enzyme	E007
 cDNA Amp Buffer	B007
 cDNA Amp Primer Mix	P003
 cDNA Amp Enzyme	E006
 TE	B009
Pouch 2 (-80°C)	
CryoCube	JW001

Variations of the Curio Seeker Dual Indexing Primer Kit: v1 and v2











To improve the workflow and ensure success for our users, we have updated the Curio Seeker Dual Indexing Primer Kit with a new version. Both v1 and v2 Dual Indexing Primer Kits support up to five 10x10 Curio Seeker tiles. Please follow the appropriate portion in this user guide for setting up the indexing PCR in section H based on the Dual Indexing Primer Kit version you are using.

Curio Seeker Dual Indexing Primer Kit v1 components

Component	Part number
Curio Seeker Dual Indexing Primer Kit v1 (-20°C)	K003
 Index Primer F1	P005
 Index Primer F2	P006
 Index Primer F3	P007
 Index Primer F4	P008
 Index Primer F5	P009
 Index Primer R1	P010
 Index Primer R2	P011
 Index Primer R3	P012
 Index Primer R4	P013
 Index Primer R5	P014

Note: Do not substitute with any other primers.

Curio Seeker Dual Indexing Primer Kit v2 components

Component	Part number
Curio Seeker Dual Indexing Primer Kit v2 (-20°C)	K006
 Index Primer F1	P016
 Index Primer F2	P017
 Index Primer F3	P018
 Index Primer F4	P019
 Index Primer F5	P020
 Index Primer R1	P021
 Index Primer R2	P022
 Index Primer R3	P023
 Index Primer R4	P024
 Index Primer R5	P025

Note: Do not substitute with any other primers.

Curio Seeker reaction chamber pouch components (4 reactions)

Component	Part number
Curio Seeker 10x10 Reaction Chamber Pouch	RC001
Curio Seeker 10x10 Reaction Chambers	TJ001
Curio Seeker Chamber Seals	TJ003

Curio Seeker 10x10 adapter

Component	Part number
Curio Seeker reaction chamber adapter	TJ002

Required reagents not included

Component	Vendor	Part number
Ethyl Alcohol	Sigma Aldrich	459844-1L
SPRIselect Reagent or sparQ PureMag beads	Beckman Coulter Quanta Bio	B23318 95196-005
Nextera XT Library Prep Kit	Illumina	FC-131-1024
Bioanalyzer High Sensitivity DNA kit or TapeStation High Sensitivity DNA D5000 ScreenTape	Agilent	5067-4626
TapeStation High Sensitivity DNA D5000 ScreenTape	Agilent	5067-5592
TapeStation High Sensitivity DNA D5000 Reagents	Agilent	5067-5593
Qubit 1X dsDNA HS Assay Kit	Thermo Fisher	Q33230

Note: Do not substitute with any other library prep kits.

Equipment

Component	Vendor	Part number
Cryostat	Leica	CM3050S
Single-channel pipette: 10 µL, 20 µL, 200 µL, and 1,000 µL	Rainin	17014388 17014392 17014391 17014382
Eight-channel or 12-channel pipette: 20 µL and 200 µL	Rainin	17013803 17013805
Mini centrifuge for 1.5 mL tubes	Major suppliers	-
Mini centrifuge for 0.2 mL tubes	Major suppliers	-
Eppendorf Centrifuge 5415 D or equivalent	Eppendorf	5425-55001
2 heat blocks for 1.5 mL tubes	Major suppliers	-
96-well PCR chiller rack or 96-Well Aluminum Block	MIDSCI Light Labs	5640-T4 A-7079
C1000 Touch Thermal Cycler	Bio-rad	1851148
Vortexer	Major suppliers	-
Tweezers	Ted Pella	58083-NM
2100 Bioanalyzer or 4200 TapeStation	Agilent Agilent	G2939BA/G2953CA G2991AA
DynaMag™-2 Magnet (1.5 mL)	Invitrogen	12321D
12-Tube Magnetic Separation Rack (PCR tube)	NEB	S1515S
Qubit 4 Fluorometer	Thermo Fisher	Q33238

Note: Although we do not expect significant differences in results between different models of thermal cycler, the protocol has been developed with C1000 Touch with a ramp rate of 3 °C/sec with a 96-well block.

Consumables


Component	Vendor	Part number
Filter pipette tips: 20 μ L, 200 μ L, and 1,000 μ L	Rainin	30389226
		30389240
		30389213
DNA LoBind 1.5 mL tubes	Eppendorf	022431021
DNA LoBind 2 mL tubes	Eppendorf	022431048
DNA LoBind 5 mL tubes	Eppendorf	0030108310
0.2 mL PCR strip tubes	USA Scientific	1402-4700
Optimal Cutting Temperature compound (OCT compound)	Sakura	4583
Qubit assay tubes	Thermo Fisher	Q32856

Tissue requirements and recommendations for assessing tissue quality

- Fresh frozen tissue
- Assess RNA quality of your tissue by collecting five 10 μ m sections and isolating RNA with the RNeasy Mini Kit (Qiagen P/N 74104) or equivalent. Analyze the RNA from your sections on an Agilent Bioanalyzer or TapeStation to derive a RNA integrity number (RIN). Good quality RNA should have a RIN value of at least 7. Although you may use samples with lower RIN scores at your own risk, it may result in lower data quality.
- Assess tissue quality by performing H&E staining on an adjacent section to the one used for the Curio Seeker workflow to provide information on tissue structural context and sectioning quality.
- Run a pilot Curio Seeker tile experiment on new tissue samples and perform shallow sequencing of 200 million reads to assess library quality.

Tips and techniques

- Use Eppendorf LoBind tubes (Eppendorf, Cat. No. 022431021 or 0030108310) in all steps where tubes are indicated.
- When using tweezers or forceps to transfer the tile, avoid direct contact with the beads. Instead, grip the glass slide from the beveled corner as shown in the example below:



- When working with multiple samples, clean tweezers or forceps between each sample by dipping the tips in 100% ethanol and wiping with Kimwipe to prevent cross contamination.
- Keep all enzymes on ice when preparing reaction mixes.
- If you are experiencing difficulties due to static charge during tissue sectioning, please follow these suggestions:
 1. *Ground yourself before sectioning.*
 2. *Change to a new pair of gloves.*
 3. *Wipe the back of the blue adhesive with 100% ethanol. Allow 15–20 seconds to dry before placing it in the cryostat.*
 4. *Regularly clean the entire cryostat chamber with 100% ethanol to reduce charge.*
 5. *Place a small piece of a dryer sheet in the corner of the cryostat.*

Protocol

A. Preparation before starting

1. Thaw the following reagents at room temperature and then keep on ice:
 - a. RT/SS Buffer
 - i. If precipitate is observed, heat the RT/SS buffer at 37°C for 5 minutes and briefly vortex before use.
 - b. dNTP
2. Place the reaction chamber adapter into the thermal cycler and set the thermal cycler to 52°C, keeping the lid completely open.



B. Tissue sectioning and hybridization to Curio Seeker tile

1. Equilibrate the fresh frozen tissue and CryoCube to -18°C in a cryostat (such as Leica CM3050S) for at least 20 minutes prior to sectioning. The optimal temperature for sectioning may vary depending on the tissue type.
2. Mount the tissue block and the CryoCube onto cutting blocks with Optimal Cutting Temperature compound (OCT compound).
3. Record the tile ID (ex: B0050_004) of the Curio Seeker tile as shown in the example on the right.



IMPORTANT: Each Curio Seeker tile is unique. The tile ID is required to retrieve the correct file for spatial barcode mapping of the sequencing data.




4. Section a 10 µm section of tissue.

5. Melt the section onto the tile in one of the two ways described below:
- OPTION 1: For precise placement of the region of interest.
 - Place the Curio Seeker tile slide in the cryostat to chill for 1 minute
 - Place the chilled Curio Seeker tile slide on the cutting stage and arrange the tissue section on top of the tile using a brush. Make sure that the region of interest is positioned directly over the tile.
 - With the tile and tissue section facing up, melt the tissue section onto the tile by moving the tile off the cryostat stage, and GENTLY placing a finger on the bottom of the slide glass, as shown in the example below. To avoid curling of the tissue, start from one side and slowly move your finger across the region rather than warming it from the center. A small brush can be used to hold the other end of the tissue flat during the initial melting from one end.



- OPTION 2: For quick placement of the region of interest.
 - Hold a **room temperature** tile in the tile holder with the tile facing down. Hover the tile over the region of interest. Keeping the tile horizontal, gently lower the tile to bring it into contact with the tissue section. The tissue section should melt onto the tile immediately.



6. Place the tile with the melted tissue section back into the cryostat.

7. Section a 30 μm section of the CryoCube. Alternatively, 60 μm can be used for tissues with higher expression levels.
8. Melt the CryoCube section onto the tile by placing a finger under the tile and moving it across the tile until the entire CryoCube section and tissue section are melted. Alternatively, briefly warm the tile with the tissue section by placing a finger underneath the tile for a few seconds, and melt the CryoCube section onto the tile with the stamping method described in step 5b.
9. Place the tile back in the Cryostat or in the tile holder on dry ice while you prepare the hybridization reaction mix in the next step, or store at -80°C for up to four days.



SAFE STOPPING POINT. Tiles can be stored at -80°C for up to 4 days in a sealed container.

10. Prepare the **hybridization reaction mix** following the table below and keep at room temperature:

Component	1 Reaction (μL)	Plus 5% overage (μL)
<input type="radio"/> Hyb buffer	380	399
<input type="radio"/> RNase inhibitor	20	21
Total	400	420

11. Write the sample name and tile ID on the side of a new reaction chamber.



12. Add 400 μL of Hybridization Reaction Mix to chamber 1.
13. Remove the tile from the cryostat, dry ice or -80°C and warm up the tile by placing a finger under the tile until the tissue has re-melted.
14. Carefully remove the tile from the blue adhesive with tweezers or forceps and place it in chamber 1 containing 400 μL of Hybridization Reaction Mix. Place the tile straight up and down in the middle of chamber 1. There should be no resistance while placing the tile in the chamber. If you encounter resistance, lightly tap the edge of the tile with the tweezer. Make sure the tile is completely submerged.
15. Seal the chambers with a chamber seal with a straight flat edge.
16. Incubate for 30 minutes at room temperature.

- Remove the remaining block of tissue and the CryoCube from the cryostat and store it at -80°C .

TIP: Cover the exposed tissue with a drop of O.C.T and freeze prior to storage to prevent desiccation of the tissue.

c. Reverse transcription

- Prepare the 1X RT Wash Buffer in a 1.5 mL tube following the table below for washing the tile prior to RT, and keep at room temperature:


Component	1 Reaction (μL)	Plus 5% overage (μL)
<input type="radio"/> RT/SS Buffer	80	84
<input type="radio"/> Nuclease-free Water	320	336
Total	400	420

- Prepare the RT Reaction Mix in a 1.5 mL tube following the table below and keep on ice:

Component	1 Reaction (μL)	Plus 5% overage (μL)
<input type="radio"/> RT/SS Buffer	80	84
<input type="radio"/> dNTP	40	42
<input type="radio"/> RNase Inhibitor	10	10.5
<input type="radio"/> RT Enzyme	20	21
<input type="radio"/> Nuclease-Free Water	250	262.5
Total	400	420

- Carefully remove the chamber seal from the reaction chamber.
- Add 400 μL of 1X RT Wash Buffer to chamber 2.
- Add 400 μL of RT Reaction Mix to chamber 3.



6. Using a pair of clean tweezers or forceps, remove the tile from the Hybridization Reaction Mix in chamber 1 and dip it in the 1X RT Wash Buffer in chamber 2 for 5 seconds.
7. Transfer the tile to chamber 3 with 400 μ L of RT Reaction Mix. Make sure the tile is completely submerged.
8. Remove the liquid from chambers 1 and 2 and discard.
9. Seal the reaction chamber with a chamber seal with a straight flat edge.
10. Incubate at room temperature for 10 minutes.
11. Place the reaction chamber onto the reaction chamber adapter in the thermal cycler that was preheated to 52°C.



12. Incubate for 30 minutes with the thermal cycler lid open.

D. Tissue clearing and Curio Seeker bead resuspension

1. If precipitate is observed, heat the TC buffer at 37°C for 5 minutes and briefly vortex before use.
2. If the room temperature in the lab is <20°C, warm the bead wash buffer to 37°C to prevent precipitation during the bead wash steps.
3. Make the Tissue Clearing Reaction Mix in a 1.5 mL tube following the table below and keep at room temperature:

Component	1 Reaction (μ L)	Plus 5% overage (μ L)
 TC Buffer	392	411.6
 TC Enzyme	8	8.4
Total	400	420

4. Take the reaction chambers out of the thermal cycler.

5. Set the thermal cycler to 37°C.
6. Carefully remove the chamber seal from the reaction chamber.
7. Add 400 µL of TC Clearing Reaction Mix into chamber 4.
8. Move the tile from chamber 3 to chamber 4.
9. Remove the liquid from chamber 3 and discard.
10. Seal the reaction chamber with a chamber seal with a straight flat edge.
11. Place the reaction chamber back onto the reaction chamber adapter in the thermal cycler that is set at 37°C.
12. Incubate for 30 minutes with the thermal cycler lid open.
13. After the incubation, carefully remove the chamber seal.
14. Transfer the tile to a 5 mL Lo-bind tube.
15. Carefully transfer the liquid from chamber 4 to the 5 mL Lo-bind tube.
16. Rinse chamber 4 with 400 µL of Bead Wash Buffer and transfer the contents to the 5 mL tube containing the tile.
17. Add an additional 500 µL of Bead Wash Buffer to the 5 mL tube containing the tile.
18. Dissociate beads from the glass slide by pipetting the wash buffer mixture directly onto the beads. Be careful not to create excess bubbles as it will make it difficult to see the tile during bead dissociation.




IMPORTANT: It is important to completely dissociate the beads from the glass slide for the region of the tile that was covered by your tissue section.

TIP: Aim the pipette tip at the bead patch when dissociating the beads from the glass slide. You should see patches of beads detach from the glass slide as you pipette. Occasionally the residual tissue and beads may still be loosely associated after detachment from the glass slide. Pipet up and down until the cloud of tissue and beads are fully dissociated.

TIP: Areas of the tile that are not covered by tissue may be harder to dissociate. For areas that are harder to dissociate, mechanically dislodge the beads from the tile by gently brushing the tip against the tile sideways while the tip is parallel to the tile and pipetting across the tile where beads remain. Incomplete dissociation of beads not covered by tissue will NOT affect data quality and performance.

19. Remove the glass from the tube with tweezers or forceps and discard.
20. Transfer the bead suspension to a new 1.5 mL tube.
21. Pellet the beads for 3 mins at 3000 x g. A white bead pellet should be visible to the eye.



22. Remove any bubbles from the top of the supernatant.
23. Carefully remove and discard the supernatant without disturbing the bead pellet as shown above.



IMPORTANT: Remove the supernatant immediately after centrifugation as the pellet may slide to the bottom of the tube after some time and become harder to visualize. When removing the supernatant, angle the tip away from the pellet and aspirate slowly to not disturb the pellet. Take care to prevent bead loss while removing the supernatant. **If the pellet begins to slide, spin for an additional 30 seconds at 3000 x g. When processing multiple samples, repeat 30 second centrifugation immediately before aspirating supernatant. .**

THIS APPLIES TO ALL BEAD WASH STEPS.

TIP: In order to retain as many beads as possible, it is not necessary to remove all the supernatant from initial washes in step. You may leave 15–20 μL . But You should leave $<10 \mu\text{L}$ of bead wash buffer before adding a reaction mix to the beads. THIS APPLIES TO ALL BEAD WASH STEPS.

24. Resuspend the bead pellet in 800 μL of Bead Wash Buffer and pellet the beads by centrifuging for 3 mins at 3000 x g.
25. *Immediately* remove and discard the supernatant. You may leave $<10 \mu\text{L}$ of Bead Wash Buffer to preserve the bead pellet.
26. Resuspend the bead pellet 400 μL of Bead Wash Buffer.








SAFE STOPPING POINT. Beads can be stored at 4°C for up to 3 days.

E. Second strand synthesis

1. Set one heat block to 95°C and another to 37°C.
2. Thaw the following reagents at room temperature and then keep on ice:

- a. RT/SS Buffer
 - b. dNTP
 - c. SS Primer
3. Keep SS Enzyme on ice.
4. Gently pipette mix the beads from the previous step 5 times.
5. Incubate the beads at 95°C for 5 minutes.
6. Prepare the Second Strand Mix in a 1.5 mL tube and keep at room temperature.

Component	1 Reaction (μL)	Plus 5% overage (μL)
 RT/SS Buffer	80	84
 dNTP	40	42
 SS Primer	4	4.2
 SS Enzyme	10	10.5
 Nuclease-free water	266	279.3
Total	400	420

7. After 5 minutes of incubation at 95°C, immediately spin beads down for 30 seconds at 3000 x g and carefully remove and discard the supernatant.







IMPORTANT: Remove supernatant immediately after centrifugation. THIS APPLIES TO ALL BEAD WASH STEPS.

TIP: If you are processing multiple samples, please leave samples at 95°C until you are ready to process them, for up to 10 additional minutes.

8. Immediately resuspend the beads in 400 μL of Second Strand Mix.
9. Incubate at 37°C for 1 hour.
10. Add 400 μL of Bead Wash Buffer.
11. Spin the beads down for 3 minutes at 3000 x g and *immediately* remove and discard the supernatant.
12. Resuspend the beads in 400 μL of Bead Wash Buffer.

F. cDNA amplification

1. Thaw the following reagents at room temperature and then keep on ice:
 - a. cDNA Amp Buffer
 - b. cDNA Amp Primer Mix
2. Keep cDNA Amp Enzyme on ice.
3. Preheat a thermocycler to 98°C and heat the lid to 105°C.
4. Prepare the cDNA Amplification Mix in a 1.5 mL tube following the table below and keep on ice:

Component	1 Reaction (μL)	Plus 5% overage (μL)
 cDNA Amp Buffer	200	210
 cDNA Amp Primer Mix	16	16.8
 cDNA Amp Enzyme	8	8.4
 Nuclease-free water	176	184.8
Total	400	420

5. Spin the beads down for 2 minutes at 3000 x *g* and *immediately* remove the supernatant.
6. Add 400 μL of the cDNA Amplification Mix to beads.
7. Split the reaction mix into 8 PCR tubes (50 μL each).
8. Pipette-mix each PCR tube before placing the tubes into the thermal cycler.

9. Immediately run the cDNA amplification program on the pre-heated Thermocycler as follows:

Temperature Ramp rate: 3°C/s	Time	Cycle
98°C	2 min	
98°C	20 sec	4 cycles (phase 1)
65°C	45 sec	
72°C	3 min	
98°C	20 sec	9* cycles (phase 2)
67°C	20 sec	
72°C	3 min	
72°C	5 min	
4°C	Hold	



IMPORTANT: If starting with a tissue section that does not cover the entire tile or has low cellular content, additional cycles for **phase 2** may be required. For example, if your tissue only covers 50% of the tile, increase the cycle number by 1-2. Add 2-3 cycles for tissues with low RNA abundance.

Fraction of tile covered by tissue	Recommended # of phase 2 cycles
>2/3 but not completely covered	10-11
2/3	11-12
1/3	12-13

G. Purification and quantification

NOTE: sparQ PureMag beads can be used for all steps where SPRI beads are indicated. Equilibrate PureMag beads to room temperature for 30 minutes prior to use.


First 0.6X bead purification

1. Prepare fresh 80% ethyl alcohol.
2. Combine the 8 reaction mixtures into one new 1.5 mL tube.
3. Measure the total volume of the combined reaction mixture and calculate the volume of SPRI beads needed: $[\text{total volume}] \times [0.6]$.
 - a. For example, if the total volume is 400 μL , you will need 240 μL of SPRI beads.
4. Vortex the SPRI beads at high speed for 30 seconds. The beads should appear homogeneous and uniform in color.
5. Add the volume of SPRI beads calculated in step 3 to the tube of combined reaction mixture.
6. Vortex to mix for 10–15 seconds.
7. Incubate at room temperature for 5 minutes.
8. Briefly centrifuge the tube.
9. Place the tube on the magnetic rack. Once the solution is clear, carefully aspirate, and discard the supernatant.
10. Keeping the tube on the magnetic stand, add 500 μL of 80% ethanol.
11. Wait 30 seconds and remove the supernatant.
12. Add 500 μL of 80% ethanol.
13. Wait 30 seconds and remove the supernatant.
14. Briefly spin the tube to collect the remaining ethanol at the bottom of the tube.
15. Place the tube back on the magnetic rack and remove the remaining ethanol.
16. Let the SPRI beads dry at room temperature until the beads appear matte (~10 minutes).
17. Remove the tube from the magnetic rack and add 50 μL of nuclease free water to the tube and pipette the beads to mix well.
18. Incubate at room temperature for 1 minute.
19. Place the tube back on the magnetic rack.
20. Once the solution is clear, transfer the supernatant to a new 0.2 mL PCR tube. Discard the used beads.

Second 0.6X bead purification

21. Vortex the SPRI beads and add 30 μ L to the tube (0.6X volume of amplification volume).
22. Vortex to mix for 10–15 seconds.
23. Incubate at room temperature for 5 minutes.
24. Briefly centrifuge the tube
25. Place the tube on the magnetic rack. Once the solution is clear, carefully aspirate, and discard the supernatant.
26. Keeping the tube on the magnetic stand, add 200 μ L of 80% ethanol.
27. Wait 30 seconds and remove the supernatant.
28. Add 200 μ L of 80% ethanol.
29. Wait 30 seconds and remove the supernatant.
30. Briefly spin the tube to collect the remaining ethanol at the bottom of the tube.
31. Place the tube back on the magnetic rack and remove the remaining ethanol.
32. Let the SPRI beads dry at room temperature until the beads appear matte (30 seconds – 2 minutes).
33. Remove the tube from the magnetic rack and add 20 μ L of nuclease free water to the tube to elute. Pipette the beads to mix well and incubate at room temperature for 1 minute.
34. Place the tube on a magnetic rack and incubate for 1 minute.
35. Transfer the supernatant to a new 0.2 mL PCR tube.
36. Quantify the cDNA products using Qubit (1X dsDNA HS assay kit) and a Bioanalyzer or TapeStation following manufacturer's guidelines (Bioanalyzer High Sensitivity DNA assay, TapeStation High Sensitivity D5000 assay or TapeStation D5000 assay). Concentrations in the range of **1 ng/ μ L** and above are acceptable. If a significant amount of primer dimer is present, you may repeat one extra round of bead purification by bringing the total volume to 50 μ L with water and following steps 21–36.

TapeStation D5000



SAFE STOPPING POINT. cDNA samples can be stored -20°C for one week before proceeding to the next step.

H. Tagmentation (with Nextera XT Library Prep Kit)

1. Preheat a thermocycler to 55°C and heat the lid to 105°C .
2. Thaw the following reagents at room temperature and then keep on ice:
 - a. Dual Indexing Primers
 - b. Nextera TD buffer
3. Keep ATM and Nextera PCR mix (NPM) on ice.
4. In a new 1.5 mL tube, add 4.8 ng of cDNA and top off to a total volume of 40 μL with nuclease-free water.
5. Add 80 μL of the Nextera TD buffer.
6. Add 40 μL of ATM.
7. Pipette to mix.
8. Split the reaction mix into 8 PCR tubes (20 μL each).
9. Briefly centrifuge the tubes.
10. Incubate at 55°C for 5 minutes.
11. After 5 minutes of incubation, *immediately* add 5 μL of Neutralization Tagment Buffer (NT) to each tube. Mix by pipette-mixing and spin down.
12. Incubate at room temperature for 5 minutes.



IMPORTANT: Identify the part number of your Curio Seeker Dual Indexing Primer Kit. If you are using the Curio Seeker Dual Indexing Primer Kit (**K003**), proceed to

Appendix A. If you are using the Curio Seeker Dual Indexing Primer Kit v2 (K006), proceed to step 13 below.

13. Add 15 µL of Nextera PCR Mix to each tube.
14. Add 5 µL of Index primer F and 5 µL of Index Primer R from the Curio Seeker Dual Indexing Primer Kit to each tube. Use the same index primers for each of the 8 partitions from the same sample.



IMPORTANT: Please ensure that each sample that will be sequenced together has a unique combination of F and R primers. One F and R primer set will be sufficient for the 8 partitions of a single tile.

15. Pipette to mix and briefly centrifuge.
16. Run the indexing PCR program according to the table below:

Temperature	Time	Cycle
72°C	3 mins	
95°C	30s	
95°C	10s	10 cycles
55°C	30s	
72°C	30s	
72°C	5 mins	
4°C	Hold	

I. Library cleanup and QC

NOTE: sparQ PureMag beads can be used for all steps where SPRI beads are indicated. Equilibrate PureMag beads to room temperature for 30 minutes prior to use.

First 0.6X bead purification

1. Prepare fresh 80% ethyl alcohol.


2. Remove the PCR tubes from the thermal cycler.
3. Combine the reaction mixture into one new 1.5 mL tube.
4. Measure the total volume of the combined reaction mixture.
5. Calculate the volume of SPRI beads needed by multiplying the total reaction mixture volume by 0.6.
 - a. For example, if the total volume is 400 μ L, you will need 240 μ L of SPRI beads.
6. Vortex the SPRIselect reagent at high speed for 30 seconds. The beads should appear homogeneous and uniform in color.
7. Add the volume of SPRI beads calculated in step 4 to the tube of amplified cDNA. Pipette Mix.
8. Vortex to mix for 10–15 seconds.
9. Incubate at room temperature for 5 minutes.
10. Briefly centrifuge the tube.
11. Place the tube on the magnetic rack. Once the solution is clear, carefully aspirate, and discard the supernatant.
12. Keeping the tube on the magnetic stand, add 500 μ L of 80% ethanol.
13. Wait 30 seconds and remove the supernatant.
14. Add 500 μ L of 80% ethanol.
15. Wait 30 seconds and remove the supernatant.
16. Briefly spin the tube to collect the remaining ethanol at the bottom of the tube.
17. Place the tube back on the magnetic rack and remove the remaining ethanol.
18. Let the SPRI beads dry at room temperature until the beads appear matte (3–5 minutes).
19. Remove the tube from the magnetic rack and add 50 μ L of TE to the tube and pipette the beads to mix well.
20. Incubate at room temperature for 1 minute.
21. Place the tubes on a magnetic rack and incubate for 1 minute.
22. Transfer the supernatant to a new 0.2 mL PCR tube.

Second 0.8X bead purification

23. Vortex the SPRIselect reagent at high speed for 30 seconds.
24. Add 40 μ L of SPRI beads to the tube (0.8X of the eluted volume).
25. Vortex to mix for 10–15 seconds.
26. Incubate at room temperature for 5 minutes.
27. Briefly centrifuge the tube.
28. Place the tube on the magnetic rack. Once the solution is clear, carefully aspirate, and discard the supernatant.
29. Keeping the tube on the magnetic stand, add 200 μ L of 80% ethanol.

30. Wait 30 seconds and remove the supernatant.
31. Add 200 μL of 80% ethanol.
32. Wait 30 seconds and remove the supernatant.
33. Briefly spin the tube to collect the remaining ethanol at the bottom of the tube.
34. Place the tube back on the magnetic rack and remove the remaining ethanol.
35. Let the SPRI beads dry at room temperature until the beads appear matte (30 seconds - 2 minutes).
36. Remove the tube from the magnetic rack and add 10 μL of nuclease free water to the tube to elute. Pipette the beads to mix well and incubate at room temperature for 1 minute.
37. Place the tube on a magnetic rack and incubate for 1 minute.
38. Transfer the supernatant to a new 0.2 mL PCR tube.
39. Quantify the cDNA products using Qubit (1X dsDNA HS assay kit) and a Bioanalyzer or TapeStation following manufacturer's guidelines (Bioanalyzer High Sensitivity DNA assay, TapeStation High Sensitivity D5000 assay or TapeStation D5000 assay). Concentrations should be $>1 \text{ ng}/\mu\text{L}$. See example library trace below:

TapeStation D5000



SAFE STOPPING POINT. Libraries can be stored -20°C before proceeding to the next step or for long term storage

J. Sequencing

- For shallow sequencing to check library quality and preview spatial expression, ~200M reads per tile is recommended.
- For deep sequencing, 1-3B reads per tile is recommended depending on the specific tissue type and tile coverage.
- Read lengths required:
 - Read 1: 50 bp
 - Index 1: 8 bp
 - Index 2: 8 bp
 - Read 2: minimum 50 bp



IMPORTANT: DO NOT PERFORM ADAPTER TRIMMING.

- Loading concentration recommendations:
 - Nextseq 1000/2000:
 - Start at 750pM final loading concentration and adjust based on sequencing quality
 - Nextseq 500/550:
 - Start at 1.8pM final loading concentration and adjust based on sequencing quality
 - Novaseq 6000:
 - Start at 250–500pM final loading concentration. Adjust based on sequencing quality
- PhiX spike-in recommendations:
 - NextSeq 1000/2000: 5% PhiX spike-in
 - NextSeq 500/550: 5% PhiX spike-in
 - NovaSeq 6000:
 - 5% PhiX spike-in when pooling with non-Curio Seeker libraries
 - 10% PhiX spike-in when pooling with only Curio Seeker libraries

- Index sequences

Index 1 primers

i7 bases for Illumina sample sheet

Index Primer F1	TAAGGCGA
Index Primer F2	CGTACTAG
Index Primer F3	AGGCAGAA
Index Primer F4	TCCTGAGC
Index Primer F5	GGACTCCT

i5 bases for Illumina Sample Sheet (NovaSeq 6000 with v1.0

Index 2 primers reagent kits, MiSeq, HiSeq 2000/2500, NextSeq 2000 (Sample Sheet v2))

Index Primer R1	TATCCTCT
Index Primer R2	AGAGTAGA
Index Primer R3	GTAAGGAG
Index Primer R4	ACTGCATA
Index Primer R5	AAGGAGTA

- Library structure



- Oligonucleotide sequences

Primers**Sequence**

SS Primer	AAGCAGTGGTATCAACGCAGAGTGANNNGGNNNB
cDNA Amp Primer Mix	CTACACGACGCTCTTCCGATCT AAGCAGTGGTATCAACGCAGAGT
Index primer F	CAAGCAGAAGACGGCATAACGAGAT-N8-GTCTCGTGGGCTCGG
Index primer R	AATGATACGGCGACCACCGAGATCTACAC-N8-ACACTCTTTCCC TACACGACGCTCTTCCGATCT

K. Informatics

FASTQ files generated by sequencing of the Curio Seeker libraries can be analyzed by the Curio Seeker bioinformatics pipeline. There are two options to access the pipeline:

- 1) To analyze your data on our cloud based analysis platform, request access at:
<https://knowledgebase.curiobioscience.com/bioinformatics/bioinfo-portal/>
- 2) To install the pipeline locally at your institution, request access at:
<https://curiobioscience.com/bioinformatics-pipeline/>.

You can download your barcode file(s) by Tile ID at:
<https://curiobioscience.com/support/barcode/>.

Troubleshooting

Issue	Likely reasons	Solution
Tissue section curling during melting step	Tissue section not flat	<p>If the tissue section is curling upwards, it may help to flip the section and flatten it with a brush to prevent further curling prior to placing the section onto the tile.</p> <p>Use a small brush to flatten the tissue section prior to melting it onto the tile. For melting, start from one edge of the tissue and slowly move across the entire tissue until the entire section has been melted onto the tile. Use a brush to hold the section flat if the section begins to curl during melting.</p>
cDNA amplification product yield is low (<1ng/ μ L)	Poor tissue or RNA quality	<p>Ensure tissue block was processed following best practices for your specific tissue type.</p> <p>Assess RNA quality of the tissue block by extracting RNA from 5 – 10 μm thick sections using the Qiagen RNeasy mini kit.</p> <p>Add additional PCR cycles for cDNA amplification if the tissue is of lower quality.</p>

	Bead loss during wash steps	<p>Ensure beads are completely dissociated from the glass slide during the tissue dissociation step.</p> <p>Check that centrifugation steps are conducted with the correct speed and length. Remove supernatant immediately after spin completes.</p> <p>Take care to not pipette any beads when removing supernatant after centrifugation steps.</p>
	Tissue section did not cover the entire tile	Additional PCR cycles may be needed if the tissue section only partially covers the tile. Start by increasing the number of cycles by 2-3 cycles and adjust based on yield and library trace.
	Tissue section too thick	Ensure tissue section is 10 µm thick
	Buffer precipitation occurred during bead wash steps after tissue clearing	Ensure the bead wash buffer and centrifuge are >20°C. If working in a lab space that is cooler, warm up the buffer to 37°C in a heat block prior to use.
	Too much bead wash buffer was left behind prior to adding the cDNA amp reaction mix to the beads	Ensure <10 µL of bead wash buffer was left behind prior to resuspending beads in the cDNA amp reaction mix.
Final library yield is too low	Input for the Nextera XT DNA prep workflow was too low	Ensure proper quantification was performed on the cDNA amplification product


Difficulty dissociating beads	Tile was not stored at the correct temperature (4°C)	Areas of the tile that are not covered by tissue may be harder to dissociate. For areas that are harder to dissociate, mechanically dislodge the beads from the tile by gently brushing the tip sideways while pipetting across the tile where beads remain. Incomplete dissociation of beads not covered by tissue will <u>not</u> affect data quality and performance.
-------------------------------	--	--

Appendix A: Indexing PCR with Curio Seeker Dual Indexing Primer Kit (**K003**)



IMPORTANT: Please confirm that you are using the Curio Seeker Dual Indexing Primer Kit (K003). The Curio Seeker Dual Indexing Primer Kit v1 (K003) has been discontinued and replaced with the The Curio Seeker Dual Indexing Primer Kit v2 (K006).

1. Ensure that the 5 minute room temperature incubation in Section *H. Tagmentation*, step 13 has been performed
2. Prepare the Nextera PCR Mix following the table below:

Component	1 Reaction (μL)
Nextera PCR mix (NPM, Illumina)	15
 Nuclease-free water	8
Total	23

3. Add 23 μL of the Nextera PCR mix to the tube.
4. Add 1 μL of Index primer F and 1 μL of Index Primer R from the Curio Seeker Dual Indexing Primer Kit to the tube.
5. Pipette to mix and briefly centrifuge.
6. Run the indexing PCR program according to the table below:

Temperature	Time	Cycle
72°C	3 mins	
95°C	30s	
95°C	10s	10 cycles
55°C	30s	
72°C	30s	
72°C	5 mins	
4°C	Hold	

7. Proceed to section *I. Library Cleanup and QC*.

For Research Use Only.

Not for use in diagnostic procedures.

Document Number CB-000578-UG

Legal Notices

© 2024 Curio Bioscience, Inc. (Curio Bioscience). All rights reserved. Display, reproduction, distribution, or creation of derivative works of all or any portion of this document without the express written consent of Curio Bioscience, is strictly forbidden.

All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR HUMAN USE OR USE IN DIAGNOSTIC PROCEDURES.

You are solely responsible for complying with all laws, rules, and regulations applicable to the products and methods described herein and for reviewing any safety data sheets and complying with all recommended safety precautions related to any product, equipment, or methods described herein. Curio Biosciences disclaims all liability for your failure to comply with any applicable law, rule, regulation, or safety precaution.

NOTHING CONTAINED HEREIN SHALL CONSTITUTE ANY WARRANTY, EXPRESS OR IMPLIED, AS TO THE PERFORMANCE OF ANY PRODUCTS DESCRIBED HEREIN. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product early access agreement accompanying the acquisition of such products. Curio Bioscience provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein, and you use these products and protocols may be subject to restrictions imposed by those third parties. The use of products and methods described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such products or early access agreement accompanying the acquisition of such products, and/or other terms that have been agreed to in writing between Curio Bioscience and user.

Permission to use Curio Bioscience's marks is NOT granted by this document. This documentation may refer to the products or services offered by third parties by their brand name or company name solely for clarity, and Curio Bioscience does not claim any rights in those third party marks or names. Further any such reference should not be interpreted as an endorsement of any third party by Curio Bioscience, nor an endorsement of Curio Bioscience by a third party.