

NadPrep EZ DNA Library Preparation Kit v2 User Manual V2.1

NadPrep EZ DNA Library Preparation Module v2
NadPrep Universal Adapter (SI) Module (for MGI)
NadPrep Universal Adapter (MDI) Module (for MGI)

Thoroughly Read This Manual Before Operation.
For Research Use Only.
Version2.1

Statement

For research use only. Not for use in diagnostic procedures.

This instruction is intended for use with the NadPrep EZ DNA Library Preparation Kit v2.

Nanodigmbio Pte. Ltd. reserves all rights.

Revision History

Version	Data Released	Change
2.1	Nov 2023	Revise description of storage temperatures for kit and components.
2.0	Jan 2022	Upgrade from the v1 of library preparation module to v2.
1.2	Aug 2021	Add method of library construction of DNA fragments with a main distribution ranging from 300 bp to 350 bp.
1.1	Feb 2021	Add compatible Unique Dual Index Adapter Modules.
1.0	Jan 2021	Original version

Content

Statement	2
Revision History	2
Announcement	4
Introduction	5
Kit Content	5
NadPrep EZ DNA Library Preparation Module v2	5
NadPrep Universal Adapter Module (for MGI)	6
Equipment and Consumable	7
Workflow	8
Protocol	9
Step 1: Sample Pretreatment	9
Step 2: Fragmentation and End Repair & A-tailing	10
Step 3: Adapter Ligation	11
Step 4: Post-ligation Cleanup	12
Step 5: Library Amplification	12
Step 6: Post-amplification Cleanup	14
Step 7: Library Evaluation	15
Step 8: Sequencing	15
Appendix 1: Grading Standard of FFPE DNA	16
Appendix 2: Size Selection	17
Appendix 3: Index Information	21

Announcement

Environment

- Temperature: 20-25°C; Humidity: 40-60%.
- The personal protective equipment should be equipped during the experiment period.

Temperature Control

- Avoid repeated freezing and thawing.
- Unless otherwise specified, the preparation of enzymatic reaction should be carried out on ice/cold metal bath.
- The fragmentase should be added at the end of reaction mix preparation to prevent its activity at room temperature.
- Make sure the thermal cycler reaches the designated temperature before loading any reaction tube in.

Avoid Cross-contamination

- Briefly centrifuge all the tubes to collect the contents to the bottom before opening the lids.
- Physically separate the laboratory space, equipment and consumables among different steps.
- Clean the lab area using 0.5% sodium hypochlorite once finishing the experiment.

Others

- Make sure the reagents are fully thawed and thoroughly mixed prior to use.
- Unless otherwise specified, all the centrifuge steps are conducted at room temperature (20-25°C).
- Unless otherwise specified, all mixing steps listed as "mix thoroughly" should be performed by either vortexing for 10 sec or pipetting up and down for 10 times, and then briefly centrifuge to collect the contents.

Introduction

NadPrep EZ DNA Library Preparation Kit v2 is designed for preparation of high-quality sequencing libraries from double-stranded DNA (dsDNA) on MGI platforms. To simplify the experimental process, multiple processes were applied in one single step, including the fragmentation, end repair and adapter ligation. This A-T ligation-based kit applies to the whole genome sequencing with DNA input ranging from 5 to 500 ng, and is compatible with hybridization capture based targeted sequencing.

The kit contains magnetic beads for DNA purification, and is optimized on size selection with NadPrep Universal Adapter (SI) Module (for MGI) and NadPrep Universal Adapter (MDI) Module (for MGI). All the components have been verified by quality control with a stable and remarkable performance.

Kit Content

NadPrep EZ DNA Library Preparation Kit v2 consists of Library Preparation Module and Adapter Module, which can be used in multiple combinations as required.

NadPrep EZ DNA Library Preparation Module v2

Catalog#	Item	Package/Storage
1002601	NadPrep EZ DNA Library Preparation Module v2, 24 rxn	Box1 / -25 ~ -15°C Box2 / 2 ~ 8°C
1002602	NadPrep EZ DNA Library Preparation Module v2, 96 rxn	Box1 / -25 ~ -15°C Box2 / 2 ~ 8°C

Package#	Component	Volume 1002601 24 rxn	Volume 1002602 96 rxn	Storage
Box1	FERA Buffer	150 µL	575 µL	-25 ~ -15°C
	FERA Enzyme	150 µL	575 µL	-25 ~ -15°C
	Enhancer Buffer	250 µL	950 µL	-25 ~ -15°C
	Ligation Buffer	750 µL	2×1,500 µL	-25 ~ -15°C
	DNA Ligase	60 µL	230 µL	-25 ~ -15°C
	2X HiFi PCR Master Mix	720 µL	2×1,600 µL	-25 ~ -15°C
	Nuclease Free Water	2.5 mL	8 mL	2 ~ 8°C
	TE Solution	1 mL	4 mL	2 ~ 8°C
Box2	NadPrep SP Beads	4 mL	15 mL	2 ~ 8°C

NadPrep Universal Adapter Module (for MGI)

Universal Adapter Module with Single Index

Catalog#	Item	M-Adapter (SI) Volume	M-Index Primer Mix (SI)#	M-Index Primer Mix (SI) Volume	Storage
1003611	NadPrep Universal Adapter (SI) Module Set A1 (for MGI), 24 rxn	60 µL	1-12	12 X 15 µL (Tube)	-25 ~ -15°C
1003621	NadPrep Universal Adapter (SI) Module Set B1 (for MGI), 96 rxn	230 µL	1-24	24 X 25 µL (Tube)	-25 ~ -15°C
1003622	NadPrep Universal Adapter (SI) Module Set B2 (for MGI), 96 rxn	230 µL	25-48	24 X 25 µL (Tube)	-25 ~ -15°C
1003620	NadPrep Universal Adapter (SI) Module Set D1 (for MGI), 480 rxn	1,100 µL	1-96	96 X 30 µL (Plate)	-25 ~ -15°C

Universal Adapter Module with Unique Dual Index

Catalog#	Item	M-Adapter (MDI) Volume	M-Index Primer Mix (MDI)#	M-Index Primer Mix (MDI) Volume	Storage
1003711	NadPrep Universal Adapter (MDI) Module Set A1 (for MGI), 24 rxn	60 µL	1-12	12 X 15 µL (Tube)	-25 ~ -15°C
1003721	NadPrep Universal Adapter (MDI) Module Set B1 (for MGI), 96 rxn	230 µL	1-24	24 X 25 µL (Tube)	-25 ~ -15°C
1003722	NadPrep Universal Adapter (MDI) Module Set B2 (for MGI), 96 rxn	230 µL	25-48	24 X 25 µL (Tube)	-25 ~ -15°C
1003725	NadPrep Universal Adapter (MDI) Module Set B3 (for MGI), 96 rxn	230 µL	49-72	24 X 25 µL (Tube)	-25 ~ -15°C
1003726	NadPrep Universal Adapter (MDI) Module Set B4 (for MGI), 96 rxn	230 µL	73-96	24 X 25 µL (Tube)	-25 ~ -15°C
1003731	NadPrep Universal Adapter (MDI) Module Set D1 (for MGI), 480 rxn	1,100 µL	1-96	96 X 30 µL (Plate)	-25 ~ -15°C
1003732	NadPrep Universal Adapter (MDI) Module Set D2 (for MGI), 480 rxn	1,100 µL	97-192	96 X 30 µL (Plate)	-25 ~ -15°C
1003733	NadPrep Universal Adapter (MDI) Module Set D3 (for MGI), 480 rxn	1,100 µL	193-288	96 X 30 µL (Plate)	-25 ~ -15°C
1003734	NadPrep Universal Adapter (MDI) Module Set D4 (for MGI), 480 rxn	1,100 µL	289-384	96 X 30 µL (Plate)	-25 ~ -15°C

⚠ Note: The library requires single strand circularization before sequencing on the MGI platform. This instruction does not contain the library circularization.

Equipment and Consumable


Equipment

Item	Description
Digital electrophoresis	Agilent 2100 Electrophoresis Bioanalyzer® system (Cat # G2939AA) Agilent 2200 TapeStation system (Cat # G2965AA) Bioptic Qsep100 capillary gel electrophoresis system or equivalent
Pipettor	General laboratory supplier
Thermal cycler	General laboratory supplier
Benchtop centrifuge	General laboratory supplier
Microcentrifuge	General laboratory supplier
Vortex mixer	General laboratory supplier
Magnet stand	Thermo Fisher DynaMag™ -2 Magnet (Cat # 12321D) Thermo Fisher DynaMag™ -96 Side Magnet (Cat # 12331D) BORTHEE-96 Side Magnet (Cat # MAG-96-11) or equivalent
Qubit fluorometer	Thermo Fisher Qubit™ 3.0 Fluorometer (Cat # Q33216) or equivalent
Timer	General laboratory supplier

Consumable

Item	Description
Absolute Ethanol	General laboratory supplier, analytical grade
0.2 mL PCR Tube	Axygen MAXYMum Recovery™ PCR Tubes, 0.2 mL flat cap (Cat # PCR-02-L-C) or equivalent

Workflow



NadPrep M-Adapter (MDI)

Safe Stopping Point

Option Optional Step

*The schematic diagram takes the unique dual index (MDI) adapter as an example.

Protocol

Sample Requirement

Sample Quality: gDNA without serious degradation (tested by gel electrophoresis analysis); FFPE DNA with one clear strip with about 15 kb in size (refer to [Appendix 1: Grading Standard of FFPE DNA](#)).

Input Amount: 5-500 ng, > 0.5 ng/μL (quantified by Qubit™ 3.0 Fluorometer).

Sample Purification: High quality DNA input, with $OD_{260}/OD_{280}=1.8-2.0$ and $OD_{260}/OD_{230}=2.0-2.5$.

! Note: EDTA inhibits the activity of fragment enzyme. If the extracted DNA sample contains EDTA, the sample must be pretreated according to [Step 1](#), otherwise the enzymatic digestion effect will be affected.

Step 1: Sample Pretreatment

Before the experiment, please confirm the sample solvent:

1. If the DNA sample contains EDTA, 1.8X NadPrep SP Beads is recommended to be used for purification of the sample, and library preparation shall be conducted after the sample is eluted with Nuclease Free Water. If purification is not carried out, the final concentration of EDTA in the system shall be calculated according to the following table, or the DNA system shall be configured by adding corresponding volumes of Enhancer Buffer and Nuclease Free Water according to the sample volume containing EDTA (1 mM);
2. Sample pretreatment is not required for gDNA or FFPE DNA samples with purity meeting the requirements and without EDTA; this step can be skipped.

Example: the sample solvent contains 1 mM EDTA, and the sample concentration is 100 ng/μL, and 250 ng input DNA is used for fragmentation, then the final concentration of EDTA in the fragmented DNA system is calculated as $250 \text{ ng}/100 \text{ (ng/μL)} \times 1 \text{ mM}/50 \text{ μL} = 0.05 \text{ mM}$; or the volume of sample containing EDTA (1 mM) is calculated as $250 \text{ ng}/100 \text{ (ng/μL)} = 2.5 \text{ μL}$.

For samples containing EDTA, the method of configuring DNA system with Enhancer Buffer is as follows:

Final Concentration of EDTA in the System	Volume of Sample Containing EDTA (1 mM)	Volume of Enhancer Buffer	Volume of Nuclease Free Water
0.015-0.05 mM	0.75-2.5 μL	0.5 μL	37-38.75 μL
0.1 mM	5 μL	1.5 μL	33.5 μL
0.2 mM	10 μL	4 μL	26 μL
0.3 mM	15 μL	5 μL	20 μL
0.4 mM	20 μL	6 μL	14 μL
0.5 mM	25 μL	7 μL	8 μL
0.6 mM	30 μL	8 μL	2 μL

! Note: When the final concentration of EDTA in the system is greater than 0.6 mM, excessive volume of Enhancer Buffer needs to be added, which may lead to the decline in data quality.

Step 2: Fragmentation and End Repair & A-tailing

1. Thaw FERA Buffer at room temperature. Mix thoroughly and centrifuge to collect the contents. Keep the tube on ice.
2. Thaw FERA Enzyme on ice. Briefly vortex and centrifuge to collect the contents.
3. Set up each reaction in a 0.2 mL PCR tube on ice as follows:

DNA System  ¹	40 µL
FERA Enzyme	5 µL
FERA Buffer  ²	5 µL
Total	50 µL

 **Note ¹:** If the input DNA is less than 40 µL, add Nuclease Free Water to a total volume of 40 µL. When the sample Input is ≤ 10 ng, add 1 µL of Enhancer Buffer to obtain the expected fragment size. At this time, the input DNA should be ≤ 39 µL.

 **Note ²:** When configuring FERA Enzyme and FERA as Mix for addition, pipette continuously for 3-5 times. Pipette slowly and carefully to ensure the expected volume.

4. Mix thoroughly and centrifuge to collect the contents.
5. Incubate the tube in the thermal cycler programmed as follows. Make sure the thermal cycler is stabilized at 25°C before loading.

25°C	See the following table for enzymatic digestion time
65°C	30 min
10°C	Hold

 **Note:** The heated lid should be set at 70°C.

Expected Insert Size	Recommended enzymatic digestion time
200 bp	30 min
250 bp	25 min
300 bp	20 min
350 bp	15 min

 **Note:** The above enzymatic digestion time is applicable to blood gDNA samples without degradation. If the sample quality varies, the enzymatic digestion size may be different. It is suggested to adjust the fragmented enzymatic digestion time according to the reference time. The fragmented enzymatic digestion time shall not be less than 10 min, as insufficient enzymatic digestion time will affect the efficiency of library preparation.

6. Store the unused reagents at -25 ~ -15°C.

Option **Optional steps:** For gDNA or FFPE DNA samples, fragments size can be selected according to downstream analysis requirements. For specific methods, see [Appendix 2 1.1](#).

Step 3: Adapter Ligation

1. Thaw Ligation Buffer at room temperature, mix thoroughly and centrifuge to collect the contents. Keep the tube on ice.
! Note: The Ligation Buffer is viscous. Please pipette slowly and carefully to ensure the expected volume.
2. Thaw DNA Ligase on ice. Briefly vortex and centrifuge to collect the contents.
3. Remove the tube in **Step 2** from the thermal cyclor and set up each reaction on ice as follows:

Products of Step 2	50 µL
NadPrep M-Adapter (SI)/NadPrep M-Adapter (DI) (as appropriate) ! 1	2 µL
Ligation Buffer ! 2	26 µL
DNA Ligase ! 2	2 µL
Total	80 µL

! Note 1: Make sure to add NadPrep M-Adapter (SI)/NadPrep M-Adapter (DI) before other reagents to avoid the adapter self-ligation.

! Note 2: If working on multiple samples, prepare Ligation Master Mix in advance with Ligation Buffer and DNA Ligase.

The dilution ratio of NadPrep M-Adapter(SI)/NadPrep M-Adapter (DI) according to the amount of input DNA is recommended as follows:

Input DNA	Adapter Concentration	Dilution Ratio	Adapter Insert Molar Ratio
500 ng	15 µM	*	10:1-6:1
250 ng	15 µM	*	20:1-12:1
100 ng	15 µM	*	50:1-30:1
50 ng	15 µM	*	100:1-60:1
25 ng	15 µM	*	200:1-100:1
10 ng	15 µM	2.5	200:1-100:1
5 ng	15 µM	5	200:1-100:1

! Note: When the input DNA ≤ 10 ng, the adapter should be diluted prior to use. For the distribution of other fragments, the dilution ratio can be adjusted on the premise of ensuring the molar ratio.

4. Mix thoroughly and centrifuge to collect the contents.
5. Incubate the tube in the thermal cyclor and run the Cycling Program III set as follows. Make sure the thermal cyclor is stabilized at 20°C before loading.

20°C	15 min
4°C	Hold

! Note: The heated lid should be left open.

6. Store the unused reagents at -25 ~ -15°C.

Step 4: Post-ligation Cleanup

Option For gDNA or FFPE DNA samples, fragments can be screened according to downstream analysis requirements. For specific methods, see [Appendix 2 1.2](#).

1. Equilibrate the NadPrep SP Beads to room temperature for a minimum of 30 min. Thoroughly vortex the beads prior to use.
2. Add 40 μ L of NadPrep SP Beads to the same tube from [Step 3](#). Mix thoroughly and then incubate the tube at 25°C for 5-10 min.
3. Briefly centrifuge and then place the tube on the magnet stand until the solution is clear (5 min). Carefully remove and discard the supernatant without disturbing the beads.

! Note: The incubation time should be adjusted according to the actual performance of various magnet stands to ensure that the solution is completely clarified.

4. Keep the tube on the magnet stand, add 150 μ L of 80% freshly prepared ethanol without disturbing the beads. Incubate the tube on the magnet stand at room temperature for 30 sec. Carefully remove and discard the ethanol.
5. Repeat step 4 once.
6. Briefly centrifuge and then place the tube on the magnet stand. Carefully remove all residual ethanol without disturbing the beads by using a 10- μ L pipette tip.
7. Open the lid, dry the beads on the magnet stand at room temperature for 5 min.

! Note: Over-dried beads may be more difficult to suspend.

8. Remove the tube off the magnet stand. Add 20 μ L of Nuclease Free Water into the tube with beads inside for [Step 5](#).

Step 5: Library Amplification

1. Thaw 2X HiFi PCR Master Mix and NadPrep M-index Primer Mix (SI)/NadPrep M-index Primer Mix (MDI) on ice. Mix thoroughly and centrifuge to collect the contents.
2. Set up each library amplification reaction in a 0.2 mL PCR tube on ice as follows:

Purified post-ligation from Step 4	20 μ L
2 X HiFi PCR Master Mix	25 μ L
NadPrep M-Index Primer Mix (SI)/NadPrep M-Index Primer (MDI)	5 μ L
Total	50 μL

3. Mix thoroughly and centrifuge to collect the contents.
4. Incubate the tube in the thermal cycler and run the Cycling Program IV set as follows, with the heated lid set at 105°C.

98°C	2 min	1 cycle
98°C	15 sec	Varies (refer to the following table)
60°C	30 sec	
72°C	30 sec	
72°C	2 min	1 cycle
4°C	Hold	-

Input DNA	Recommended PCR Cycles of NadPrep M-Adapter (SI)/ NadPrep M-Adapter (DI)	
	Library Yield $\geq$ 500 ng	Library Yield $\geq$ 1,000 ng
500 ng	3	3
250 ng	3	4
100 ng	4-5	5-6
50 ng	5-6	6-7
25 ng	6-7	7-8
10 ng	8-9	9-10
5 ng	9-10	10-11

M Important notes:

The selection of PCR cycles shall take into account the following:

- If double-sided size selection has been performed, one or two additional cycles are necessary.
- The recommended cycles of DNA samples from different grades of FFPE is shown in Appendix 1: Grading Standard of FFPE Samples;
- If the library is constructed for hybridization capture-based targeted sequencing, the library yield should be twice as much as the capture input by adjusting the PCR cycles;
- To obtain 500 ng of library yield, the DNA input should be less than 250 ng.

Step 6: Post-amplification Cleanup

1. Add the corresponding volume of NadPrep SP Beads into the Library Amplification tube. Mix thoroughly and then incubate the tube at 25°C for 5-10 min.

Primer Mix Applied to Step 5	Volume of NadPrep SP Beads
NadPrep M-Index Primer Mix (SI)	60 µL
NadPrep M-Index Primer Mix (MDI)	50 µL

2. Briefly centrifuge and then place the tube on a magnet stand until the solution is clear (~5 min). Carefully remove and discard the supernatant without disturbing the beads.

! Note: The incubation time should be adjusted according to the actual performance of various magnet stands to ensure that the solution is completely clarified.

3. Keeping the tube on the magnet stand, add 150 µL of freshly prepared 80% ethanol without disturbing the beads. Incubate the tube on the magnet stand at room temperature for 30 sec. Carefully remove and discard the ethanol.

4. Repeat step 3 once.

5. Briefly centrifuge and then place the tube on the magnet stand. Carefully remove all residual ethanol without disturbing the beads by using a 10-µL pipette tip.

6. Dry the beads on the magnet stand at room temperature for 5 min.

! Note: Over-dried beads may be more difficult to suspend.

7. Remove the tube off the magnet stand. Thoroughly resuspend the beads in 20 µL of TE solution. Incubate the tube at 25°C for 5 min.

8. Briefly centrifuge and place the tube on the magnet stand. Incubate for 2 min until the solution is clear. Transfer the clear supernatant to a new 0.2 mL PCR tube. Avoid transfer of NadPrep SP Beads.



Safe Stopping Point: The DNA library can be stored at 2 ~ 8°C for one week or -25 ~ -15°C for longer.

Step 7: Library Evaluation

1. The Qubit™ fluorometer or qPCR is recommended for library quantification.
2. The Bioanalyzer® 2100 (Agilent) or Qsep 100 (Bioptic) is recommended for size distribution analysis.
3. The quality control standards include:
 - a. No adapter dimers or free adapters.
 - b. Relatively concentrated size distribution.
 - c. Suitable for MGI platform in size.




Fig 1. Example of size distribution of library. Use NadPrep EZ DNA Library Preparation Module V2 coupled with **A.** NadPrep Universal Adapter (SI) Module (for MGI) and **B.** NadPrep Universal Adapter (MDI) Module (for MGI) for size distribution of library (Agilent 2100, DNA 1000 Kit) after the library preparation of gDNA samples (with average insert fragment of 250 bp), respectively.

Step 8: Sequencing

Cyclize the DNA library with the circularization steps before loading on MGI platforms.

For targeted sequencing, please refer to the hybridization capture protocol (for MGI) from Nanodigmbio (www.nanodigmbio.com) and then perform the circularization steps.

Appendix 1: Grading Standard of FFPE DNA

FFPE DNA was classified into three grades. The corresponding experimental parameters were recommended as follows:

Grade	Grading Standard by Gel Electrophoresis	Recommended Input	Recommended Fragmentation Time (250 bp)	Recommended PCR Cycle (250 bp)
FFPE B ⁺	One clear and slight trailing strip with about 15 kb in size	≥ 100 ng	Consistent with gDNA	Consistent with gDNA
FFPE B	One indistinct strip with about 15 kb in size, with medium diffusion	≥ 100 ng	Reduce five minutes or more than gDNA	One or two extra cycles more than gDNA
FFPE C	Multiple strips ranging from 200 bp to 2,500 bp, with severe diffusion	≥ 100 ng	Reduce ten minutes or more than gDNA	Two or three extra cycles more than gDNA
FFPE D	Stripes mainly ranging from 250 bp to 1,000 bp, with severe diffusion	≥ 100 ng	Reduce fifteen minutes or more than gDNA	Three or four extra cycles more than gDNA




Fig 2. Example of size distribution of library. Use NadPrep EZ DNA Library Preparation Module v2 coupled with NadPrep Universal Adapter (MDI) Module (for MGI) for size distribution of library (Agilent 2100, DNA 1000 Kit) after the library preparation of FFPE samples (with average insert fragment of 250 bp).

Appendix 2: Size Selection

The NadPrep DNA Universal Library Preparation Kit v2 is compatible with a variety of size selection modes. Please choose the appropriate size selection modes according to the sample type, application direction, workflow and sequencing requirements.

Size Selection	Applicable Situation	Advantage	Disadvantage
Step 2 Option Double-sided size selection after fragmentation	Sufficient DNA input; Broad size distribution of DNA fragments after fragmentation	Narrow size distribution, high ligation efficiency and uniform sequencing data	Loss of sample
Step 4 Option Double-sided size selection after adapter ligation	Sufficient DNA input; Broad size distribution of libraries after adapter ligation	Narrow size distribution, Uniform sequencing data	Strict double-sized selection parameters based on library size

! Note: Double-sided size selection during library construction can generate even size distribution and is therefore beneficial to data evenness after multiplexed capture and/or multiplexed sequencing. However, loss of sample is unavoidable due to extra handling steps. Please refer to the table to choose appropriate method.

1.1 Double-sided Size Selection after DNA Fragmentation (Optional Step after Step 2)

The volume of beads for double-sided size selection of fragmented product of **Step 2** are recommended as follows:

M Important note:

- The average size of the fragmented DNA should coincide the predicted size, to avoid the aggravation of sample loss.

Expected DNA Size Distribution ¹ !	Sample Volume	Volume of NadPrep SP Beads (V1)	Volume of NadPrep SP Beads (V2) ² !
200 bp	50 µL	45 µL	50 µL
250 bp	50 µL	40 µL	50 µL
300 bp	50 µL	35 µL	50 µL
350 bp	50 µL	30 µL	50 µL

! Note ¹: For WGS, it is recommended to select DNA fragments with a mean length of 300-350 bp. For targeted sequencing, it is recommended to select DNA fragments with a mean length of 200-250 bp.

! Note ²: The beads in V2 should be concentrated. Please refer to **Appendix 1.1** step 3 as follows.

- Equilibrate the NadPrep SP Beads to room temperature for a minimum of 30 min. Thoroughly vortex the beads for 45 sec at high-speed.
- Transfer 50 µL of fragmented DNA in **Step 2 Fragmentation and End Repair & A-tailing** into a new 0.2 mL PCR tube. Add V1 volume of NadPrep SP Beads into the tube. Mix thoroughly and then incubate at 25°C for 5-10 min.

! Note: If the fragmented DNA is less than 50 µL, add Nuclease Free Water to a total of 50 µL.

- Prepare V2 volume of NadPrep SP Beads in a new 0.2 mL PCR tube. Briefly centrifuge and then place the tube on a magnet stand until the solution is clear (~5 min). Carefully remove and discard 40 µL of the supernatant without disturbing the beads. Remove the PCR tube off the magnet stand and mix the residues thoroughly.

❗ Note: The incubation time should be adjusted according to the actual performance of various magnet stands to ensure that the solution is completely clarified.

4. Briefly centrifuge and then place the PCR tube in [Appendix 1.1](#) step 2. on the magnet stand until the solution is clear (~5 min). Transfer all the supernatant (50 μ L+V1) into the PCR tube in [Appendix 1.1](#) step 3. Mix thoroughly and then incubate at 25°C for 5-10 min.

❗ Note: Keep the supernatant (50 μ L+V1) and discard the beads.

5. Briefly centrifuge and then place the tube on the magnet stand until the solution is clear (~5 min). Carefully remove and discard the supernatant without disturbing the beads.

6. Keeping the tube on the magnet stand, add 150 μ L of freshly prepared 80% ethanol without disturbing the beads. Incubate the tube on the magnet stand at room temperature for 30 sec. Carefully remove and discard the ethanol solution.

7. Repeat step 6 once.

8. Briefly centrifuge and then place the tube on the magnet stand. Carefully remove all residual ethanol without disturbing the beads by using a 10- μ L pipette tip.

9. Dry the beads on the magnet stand at room temperature for 5 min.

❗ Note: Over-dried beads may be more difficult to suspend.

10. Remove the tube off the magnet stand. Thoroughly resuspend the beads in 43 μ L of Nuclease Free Water. Incubate the tube at 25°C for 2 min.

11. Briefly centrifuge and then place the tube on the magnet stand. Incubate for 2 min until the supernatant is clear (collect 1 μ L of supernatant for quantification). Transfer the supernatant (~40 μ L) containing DNA to a new 0.2 mL PCR tube. Keep the tube on ice and directly perform [Step 3 Adapter Ligation](#).




Fig 3. DNA size distribution with double-sided size selection after fragmentation (Covaris™) by using Agilent 2100 Analyzer.

1.2 Double-sided Size Selection after Adapter Ligation (Optional Step for Step 4)

The volume of beads for double-sided size selection of ligated product of **Step 3** are recommended as follows:

M Important note:

- This step is not recommended for input DNA ≤ 50 ng or DNA with low quality.

Expected DNA Size Distribution	Volume of NadPrep SP Beads (V1)	Volume of NadPrep SP Beads (V2)
200 bp	34 μ L	32 μ L
250 bp	26 μ L	30 μ L
300 bp	22 μ L	28 μ L
350 bp	20 μ L	25 μ L

1. Equilibrate the NadPrep SP Beads to room temperature for a minimum of 30 min. Thoroughly vortex the beads for 45 sec at high-speed.
2. Add V1 volume of NadPrep SP beads into the PCR tube of **Step 3 Adapter Ligation**. Mix thoroughly and then incubate at 25°C for 5-10 min.
3. Briefly centrifuge and then place the tube on the magnet stand until the solution is clear (~5 min). Carefully remove and discard the supernatant without disturbing the beads.

❗ Note: The incubation time should be adjusted according to the actual performance of various magnet stands to ensure that the solution is completely clarified.

4. Keeping the tube on the magnet stand, add 150 μ L of freshly prepared 80% ethanol solution without disturbing the beads. Incubate the tube on the magnet stand at room temperature for 30 sec. Carefully remove and discard the ethanol solution.
5. Repeat step 4 once.
6. Briefly centrifuge and then place the tube on the magnet stand. Carefully remove all residual ethanol solution without disturbing the beads by using a 10- μ L pipette tip.
7. Dry the beads on the magnet stand at room temperature for 5 min.

❗ Note: Over-dried beads may be more difficult to suspend.

8. Remove the tube off the magnet stand. Thoroughly resuspend the beads in 50 μ L of Nuclease Free Water. Incubate the tube at 25°C for 2 min.
9. Add V2 volume of beads and then mix thoroughly. Incubate at 25°C for 5-10 min.
10. Briefly centrifuge and then place the tube on the magnet stand until the solution is clear (~5 min). Transfer all the supernatant (50 μ L+V2) to a new 0.2 mL PCR tube.

❗ Note: Keep the supernatant (50 μ L+V2) and discard the beads.

11. Add 30 μ L of NadPrep SP Beads to the supernatant (from Step 10). Mix thoroughly and incubate at 25°C for 5-10 min.
12. Briefly centrifuge and then place the tube on the magnet stand until the solution is clear (~5 min). Carefully remove and discard the supernatant without disturbing the beads.
13. Keeping the tube on the magnet stand, add 150 μ L of freshly prepared 80% ethanol solution without disturbing the beads. Incubate the tube on the magnet stand at room temperature for 30 sec. Carefully remove and discard the ethanol solution.
14. Repeat step 13 once.

15. Briefly centrifuge and then place the tube on the magnet stand. Carefully remove all residual ethanol without disturbing the beads by using a 10- μ L pipette tip.

16. Dry the beads on the magnet stand at room temperature for 5 min.

! Note: Over-dried beads may be more difficult to suspend.

17. Remove the tube off the magnet stand. Add 20 μ L of Nuclease Free Water into the tube with beads inside. Keep the tube on ice and directly perform **Step 5 Library amplification**.




Fig 4. Library size distribution with double-sided size selection after ligation by using Agilent 2100 Analyzer.

Appendix 3: Index Information

The index sequences of NadPrep M-Index Primer Mix in NadPrep Universal Adapter (SI) Module (for MGI) are listed as follows:

Index# (96-well Plate Position)	Index Sequence	Index# (96-well Plate Position)	Index Sequence
1 (A1)	TAGGTCCGAT	41 (A6)	TTAGATGCAT
2 (B1)	GGACGGAATC	42 (B6)	GTCCAGAGCT
3 (C1)	CTTACTGCCG	43 (C6)	CACGTGATAG
4 (D1)	ACCTAATTGA	44 (D6)	CCACTAGTCC
5 (E1)	TTCGTATCCG	45 (E6)	TGGACTTGGC
6 (F1)	GGTAACGAGC	46 (F6)	GCTTGACAGG
7 (G1)	CAACGTATAA	47 (G6)	AAGACCTCTA
8 (H1)	ACGTCGCGTT	48 (H6)	AGTTGCCATA
9 (A2)	TTCTGCTAGC	49 (A7)	ATGTACGCAG
10 (B2)	AGGAAGATAG	50 (B7)	TTAATGAGAT
11 (C2)	GCTCTTGCTT	51 (C7)	TGCGCCACTT
12 (D2)	CAAGCACGCA	52 (D7)	CATTAAGGCC
13 (E2)	CGGCAATCCG	53 (E7)	CCGCCTCAGA
14 (F2)	ATCAGGATTC	54 (F7)	AATCGGCTCG
15 (G2)	TCATTCCAGA	55 (G7)	GCCGGTTATC
16 (H2)	GATGCTGGAT	56 (H7)	GGAATATTGA
17 (A3)	GTGAGTGATG	57 (A8)	ATTCAACGGA
18 (B3)	GAGTCAGCTG	58 (B8)	AACTGTACTG
19 (C3)	TGTCTGCGAA	59 (C8)	GTACCTCAAT
20 (D3)	ATTGGTACAA	60 (D8)	GACTTCTAAT
21 (E3)	CGATTGTGGT	61 (E8)	TGAAGCGTTG
22 (F3)	ACAGACTTCC	62 (F8)	CGTGCGATCC
23 (G3)	TCCACACTCT	63 (G8)	TCGGAAGGCA
24 (H3)	CACCACAAGC	64 (H8)	CCGATGTCGC
25 (A4)	TAGAGGACAA	65 (A9)	ACTTAGAATG
26 (B4)	CCTAGCGAAT	66 (B9)	TCCAAGCCTG
27 (C4)	GTAGTCATCG	67 (C9)	AGACGATGAT
28 (D4)	GCTGAGCTGT	68 (D9)	CTCACAAGAC
29 (E4)	AACCTAGATA	69 (E9)	CGTTCCTACT
30 (F4)	TTGCCATCTC	70 (F9)	GTGGTTGTGA
31 (G4)	AGATCTTGCG	71 (G9)	GAAGGCCTGC
32 (H4)	CGCTATCGGC	72 (H9)	TAGCTTGCCA
33 (A5)	GCAACGATGG	73 (A10)	GACAATGCTC
34 (B5)	TAATCGTTCA	74 (B10)	GCTAATCACA
35 (C5)	GTTTCGCTCTA	75 (C10)	AGTCCATAGG
36 (D5)	TCTCACACAT	76 (D10)	CTATCGCCTA
37 (E5)	CTGTTAGGAT	77 (E10)	ATCGTGGTCT
38 (F5)	CGCAGACGCG	78 (F10)	TGGCTAATAC
39 (G5)	AAGGATCATC	79 (G10)	CAGTGCAGAG
40 (H5)	AGCGTTGAGC	80 (H10)	TCAGGCTGGT

Index# (96-well Plate Position)	Index Sequence	Index# (96-well Plate Position)	Index Sequence
81 (A11)	ATACTCACGC	89 (A12)	CACCATGTCT
82 (B11)	ATGCTCCGCG	90 (B12)	ATATGTCTGG
83 (C11)	TGTGAACTTG	91 (C12)	AAGGAAGCGT
84 (D11)	GAGAGGTGCT	92 (D12)	TCAAGACGTC
85 (E11)	TGCACTGTAA	93 (E12)	CCGCTCAGTA
86 (F11)	GCCTAGGCAA	94 (F12)	GGTGTGTACA
87 (G11)	CCATCATAGC	95 (G12)	TTCACGTAAG
88 (H11)	CATGGTAATT	96 (H12)	GGTTCCACAC

⚠ Note: The index sequences are as same as the official index sequences on MGI platform.

The single-index sequences on MGI platform are designed to be balanced within each group. To achieve the best sequencing quality, it is recommended to use continuous indexes as indicated within one block in the table and involve at least 8 indexes for multiplexed sequencing.

The index sequences of NadPrep M-Index Primer Mix in NadPrep Universal Adapter (MDI) Module (for MGI) are listed as follows:

Index# (MDI)	Primer 1 Index	Primer 2 Index	Index# (MDI)	Primer 1 Index	Primer 2 Index
1	CGTCGATGAC	TAACACGACG	25	AGGTGCCTTA	GTGCTGTCAC
2	ATATAAGGCG	TGTTCTCTTC	26	AAGAACCAAG	TCCGAGATGC
3	GATCGTGCTC	GAGTTCACAA	27	CATACATGAC	GAAGCTCTTA
4	CAGTCTTCGG	CTGATGTCCT	28	TGCCTGGTGA	TTGAGTAAGG
5	AGAACGATCT	AGACAGTGGC	29	TCTCGGAGTT	CAACACGAAT
6	TTGGTGCATT	CTCACACTTA	30	GTCGTAGACT	ACCTCCGCTA
7	GCCGTCATAA	GCCGGTAAGT	31	GCATATACCG	AGTATATGCG
8	TCCAACCAGA	ACTGGAGGAG	32	CTAGCTTCGC	CGTTGACGCT
9	GATAGCAAGA	CAACAGTAAC	33	GACGTATCAA	CTAGCAGTTA
10	ACCGTGCTTC	ATAACGCTCA	34	CCTGCTAGGA	GCATTCTCGC
11	GCAGATGTAA	GATTGCGCCT	35	CAACTTGGCG	GATCCAATAC
12	TGTTGGAGCG	CGGTGTTGGA	36	ATGCGACCTC	TTCGATGACG
13	TTGTATCCAC	TTCCAACATC	37	GTTAAGGACG	CCGCATACTA
14	CGCACAGATG	TGCGTCATGG	38	TCCAAGTTGT	TGGAGCCGAT
15	CAACTCTCGT	GCTATAGCTT	39	AGGTCCATTC	AGCTTGCGCT
16	ATGCCATGCT	ACGGCTAGAG	40	TGATGCCAAT	AATAGGTAGG
17	AGGCAGCTTA	GTAGAGATTG	41	CCTAAGAGTT	TATGCTGAGA
18	TAGCCTAGCG	TCGAGCGTTA	42	GATACTAGCT	CGGTCTCTGT
19	ATCACGTGCG	GATGAAGGAT	43	AAGGCATCTC	GTGTACGTCG
20	CGTTATGCGC	ACGAGGTACT	44	TGGCATCTGG	GTACTAAGAG
21	CAAGGATCGA	CTTCCTACAA	45	ACACTGGAGC	AGTGGATATC
22	GTTGTCGTAT	CAATCTCGCG	46	GTATGCCACA	CAACTGCGCA
23	GCAATCAATC	TGCCTATCGC	47	CTCTTAGTAG	ACCAAGACTC
24	TCCTGACAAT	AGCTTCCAGC	48	TGCGGCTCAA	TCCAGCTCAT
49	TCACATTGCT	GATAGTAACG	65	ATCGTGGATG	GACGCGGTAT
50	AATGGCGCTC	TGAGTGGCTA	66	TGGAGATCGA	ACGAGACGTC
51	GTCTCAATGA	CCGTCATTAC	67	CCTCACAGAT	CTACTCAAGA
52	CGGATGCAAG	ATCCACCGGT	68	GAATCTCTCC	TGTTATTCCG
53	AAGCCTATTG	GCTTGTTTCA	69	GCAGACTGAC	CCGTCACTGA
54	CGCTACTGCA	AACAAGCACT	70	CTCATTAACG	TGACGCAACT
55	TCAAGAGCAT	TTGCCAGTGA	71	TGGTGAGCTT	GTTGTTGCTC
56	GTTGTGCAGC	CGAGTCAGTC	72	AATCCGCTGA	AACAAGTGAG
57	AGACAGGAAT	TCGCGATGTC	73	TCGCATCAAC	TGGTTGGAAT
58	CCTTGCCGTA	AGTGACACCA	74	AGAACAGTGA	GAACGTTCCG
59	GTCATTACGG	GACATTCAAG	75	CATGTCTCCT	ATCGCACGCA
60	TAGGCATTCC	CTATCGGTGT	76	GTCTGGAGTG	CCTAACATTC
61	AGCCAGTAGG	CCGATGACGT	77	GAGGTCTGTG	CGTGTGGATG
62	GTAAGTGAC	TTATCTCGAG	78	CTATAGACGT	TTACGTCCGG
63	TAGTCACGTT	AGCCGATACC	79	TGCAGTGACC	ACGAACATAC
64	CCTGTCACCA	GATGACGTTA	80	ACTCCACTAA	GACTCATCCA

81	GCGAAGTAGG	TGCCTGACTC	121	TCCGACAGTT	CCGTAGTTGG
82	TGCCTAACCT	CATTATCGCT	122	AGTAGTTACC	TTAGCTCGTT
83	AATGGTCTAC	GTAAGCGAGG	123	GAATCGGTGG	GATCGCACAA
84	CTATCCGGTA	ACGGCATTAA	124	CTGCTACCAA	AGCATAGACC
85	TACGCTTCAG	GACACTGCCA	125	GAGGCTGAAT	TCTAATCGCA
86	CGGAGCATCT	CCGTAGCATG	126	AGTAAGAGTC	AAGCTAACAC
87	GTA TAGATC	ATTCTATGGC	127	TCCTGACTGA	CGAGCGTAGT
88	ACTTAGCGGA	TGAGGCATAT	128	CTACTCTCCG	GTCTGCGTTG
89	ATCACTCCAT	ATACGCGCCA	129	CGAAGTACCG	GCGGATATTC
90	GATCGCAGTG	TCTTAGTGAC	130	ACCGAATAGA	ATACCGCAAG
91	CCGGAATTCC	GAGGTACATT	131	TAGTCGCTAT	TGTTGCGCGT
92	TGATTGGAGA	CGCACTATGG	132	GTTCTCGGTC	CACATATGCA
93	CACAAGGTCG	TGGCTCGCTT	133	AAGAGACACA	GCAACCTGAG
94	TCTCGCAGGA	GTCACCTAAG	134	CCACTTGGAT	TAGTAGCACT
95	GTGGTATCAT	CAATAGAGGC	135	GTCTCCTTGC	AGCCTTATTC
96	AGATCTCATC	ACTGGACTCA	136	TGTGAGACTG	CTTGGAGCGA
97	CCGTTGTGAG	GCTCTGTGTC	137	AGGATAGCCG	AGACGTGTTA
98	AATGGTGAGT	TGGTGCGCAA	138	TCTGCTATAC	CCGGTGAGAT
99	TGCAACCTTC	AACACACTGT	139	GTCTACTGGT	GTCTCACAGC
100	GTACCAACCA	CTAGATAACG	140	CAACGGCATA	TATAACTCCG
101	CCGATTTCGAT	TAGTTGGTAC	141	GTACAGCGGA	GAAGTCCTTG
102	ATTGCAGCCA	CGAGACAAC	142	CGCATAGCAC	TCCAATGCCA
103	GAATACTTGC	GTTCCACCGA	143	ACTGGTTATG	ATGTCGAGAC
104	TGCCGGAATG	ACCA GTTGTG	144	TAGTCCATCT	CGTCGATAGT
105	ATGCGGATCC	CAAGCAGCAG	145	ACCATACTCC	ACATGCGTAA
106	CCTTCATGAA	GCGAACTGGT	146	CGGTGTACAA	TGCAAGTAGC
107	GACGATCATG	TGCTGTCTTA	147	GTTGCGGAGT	CTTCCAAGCT
108	TGAATCGCGT	ATTCTGAACC	148	TAACTACTGTG	GAGGTTCTCT
109	TTCGTTTACA	CAACACCGAG	149	CAAGTAACGG	CACGTTGCAG
110	GAGACGGTTG	GTTGCAGATC	150	TCGAGTGTAT	GTGAACAAGA
111	ACACAATGGT	TGGTGGTTGT	151	GTTCCGTACC	TGTTGGTGTG
112	CGTTGCACAC	ACCATTACCA	152	AGCTACCGTA	ACACCACTCT
113	AGAGTG CATG	GTCCGATTAC	153	TCTAACGTGT	CGAGTTGAAT
114	GTTTCGCATAC	CAGTCGGATT	154	AGGCCTCCAA	GACACCTTCC
115	TCGAATTCGT	TGTGATCGCG	155	GAATTGAACC	TTGCAACGTA
116	CACTCAGGCA	ACAATCACGA	156	CTCGGATGTG	ACTTGGAACG
117	AGAAGTAGGA	TCGTGGAATT	157	GTTCTAACTC	GTCATCGCGT
118	CTTCCAGTTG	AATACAGCGC	158	CCA ACTCGCT	TAGGAATGCA
119	GAGGTGCAAC	GTCCTCTTAA	159	AGGTAGTAAG	CCTTGCTCTAC
120	TCCTACTCCT	CGAGATTGCG	160	TACGGCGTGA	AGACCGAATG

161	CCATTGAATG	CGTAGATCGC	201	GCGAATCAGC	ACTGGATCCG
162	TTCAATCGAC	TAAGCGCGTG	202	AACTTCGGTA	GAGAACAGAT
163	AGGCGCTTCT	GTGCTTGAT	203	TGTGGAATAG	TTACCTCAGA
164	GATGCAGCGA	ACCTACAACA	204	CTACCGTCCT	CGCTTGTTTC
165	CGAGAATTCA	CCAGATGTGC	205	CCATAAGAGG	TCACCTCTCG
166	GCGTTCGCTT	GATTCGTGAG	206	GTCGGTAGTC	ATCAGCGCTC
167	TTCCGTAAGC	TTGCTCACCA	207	TGTCTGCCAT	CGGTAATGAA
168	AATACGCGAG	AGCAGACATT	208	AAGACCTTCA	GATGTGAAGT
169	CGGTAATATC	TGGATCGGTC	209	TACACTCGGT	ACCAATCGTT
170	ACAGCCGGAA	ATATAGCCGG	210	ATACAGGACG	TGTTTCAGCAA
171	GTTCTGCTCG	CCTCCATTAT	211	CGTGTCTCAC	CAGGTGTAGC
172	TACAGTACGT	GACGGTAACA	212	GCGTGAATTA	GTACGCATCG
173	CGTGGAGTTG	CAATTCCTCA	213	GAACACTATC	AGTTGTCTGC
174	GTCTTCTCGA	GCGAGAAGAG	214	TCTGCGGTAA	CTAGAGGACA
175	TCACCTCAAT	TTCCAGTAGT	215	AGGTTAAGCG	GCGACCTCAT
176	AAGAAGAGCC	AGTGCTGCTC	216	CTCAGTCCGT	TACCTAAGTG
177	TCTGATGGAA	ACTGGTGCAA	217	CGCAATGAAC	GTGTTATCTC
178	AGAACGTCCG	CGCCAGAATC	218	TTGGCACTTA	TCAGACGTCTG
179	GTCTTCCATT	GTAATACTCG	219	GCATTACAGC	AACCGGAAGT
180	CAGCGAATGC	TAGTCTGGT	220	AATCGGTCGT	CGTACTCGAA
181	ACAAGCTCTA	ACTGGCTAAC	221	CCTTGGTGCA	TGTGACTACA
182	CTGCAGGAAG	CTACTGATGT	222	GAGCATCCAT	ACACTTGCGG
183	TGTTCTAGCT	GACTATGGCA	223	AGAGTCGATC	CTGTGGCTAT
184	GACGTACTGC	TGGACACCTG	224	TTCACAAATGG	GACACAAGTC
185	GTAAGGTGTT	TCTCACCAGG	225	GTATCCTCAG	TCATACTGTC
186	AAGGCTAAGC	AACGTAGGCT	226	TACCGAGTGT	GATACGGAGG
187	CGCTTCCTAA	CGATCTTCTC	227	ACTGTGCACC	CGGCTAACCA
188	TCTCAAGCCG	GTGAGGATAA	228	CGGAATAGTA	ATCGGTCTAT
189	TAGGACCGTC	AGTGGATCAA	229	GTGTCTACAG	GCCTCTGTGT
190	AGTCCTGTAA	TACATCAACG	230	TACGTATACC	TATAGGACAC
191	CTAAGGAAGG	CTACCGCGTT	231	AGTAAGCGGT	AGAGAACTGA
192	GCCTTATCCT	GCGTATGTGC	232	CCACGCGTTA	CTGCTTGACG
193	ACTCGCGGTA	ACAATCCGGC	233	TAAGCTACCA	AACTGACCGT
194	CACGAACACT	TTGTGGTATG	234	AGCTGAGGTG	GTACTTATCC
195	GTAACATATGG	CGCCATGTAA	235	CTTCTCCAGC	CGTACGTAAG
196	TGGTTGTCAC	GATGCAACCT	236	GCGAAGTTAT	TCGGACGGTA
197	CGGTGAAGTC	CTACGCGAAG	237	GTGAATTGAC	AGAAGATGCA
198	GTAACATGCAT	AAGGTGATTC	238	TACCTAGACT	GTTCCGCTTC
199	TCTCTGTTGA	GCTTATTGCT	239	AGATGCCTGG	CACTATGAGT
200	AACGACCACG	TGCACACCGA	240	CCTGCGACTA	TCGGTCACAG

241	ATTGAGCACT	CCTATAATGC	281	CGCCTGTTGT	AGTTAGTCAC
242	TGAATTGTGC	AGCTACCAAT	282	GAGAATAGCA	CTCGTTGTCTG
243	GACCGATGTG	TTACGTGCCG	283	TTAGCACATC	GCGAGCCATA
244	CCGTCCACAA	GAGGCGTGTA	284	ACTTGCGCAG	TAACCAAGGT
245	ATAGGCGTTA	TCAAGTCGCC	285	CGTGATTGGT	TAGCGGAGTG
246	TGCCATACGG	GTGTACACTA	286	ATCATGCTAC	AGCAACTAGA
247	CATTCGTGCC	CGCGTAGTGT	287	GCATCCGATG	CTAGTTCTCT
248	GCGATACAAT	AATCCGTAAG	288	TAGCGAACCA	GCTTCAGCAC
249	TTAACGCGAC	TGCAGAATCC	289	GATCCTCTCG	CATAATCCGG
250	CATGTTGTCA	ACTTAGGAGT	290	TTGAAGGAAC	GTCCGCTGAT
251	AGGCGATCTG	CTGGTCTCAA	291	CCATTCTGTA	AGATCAGTCC
252	GCCTACAAGT	GAACCTCGTG	292	AGCGGAACGT	TCGGTGAATA
253	GAGGTTAGTA	TGGCTACATA	293	TTCCACCTGG	CAAGCCTGGT
254	AGTTCCGCCT	CCATGTAGAG	294	GCATCTTCTT	GTCCTTGCCA
255	CCACAGCAAG	GTTGACGTCC	295	AATGGAGAAC	TCTAGGATTG
256	TTCAGATTGC	AACACGTCGT	296	CGGATGAGCA	AGGTAACAAC
257	CGGATACACT	TAACCACCTC	297	CATGCCTATG	CTGGAGGTCT
258	GACGATAGGA	ACGTGCTTAT	298	GCGTGGACAA	AGCTCTACAA
259	TCTCGGTCTG	GTCGTGGAGA	299	AGCAATGGCT	GATATATGGC
260	ATATCCGTAC	CGTAATAGCG	300	TTACTACTGC	TCACGCCATG
261	AGATAGCCTT	AACCGTGGTA	301	AACGCCGTCA	TACCGCACAT
262	CTTACCATGC	GTTGTACTAC	302	GCGCTGAATC	GCAGTACACG
263	TCGGTTGGAG	TCATCGAAGG	303	CTAAGTTGGT	CGTACGGTTA
264	GACCGATACA	CGGAACTCCT	304	TGTTAACCAG	ATGTATTGGC
265	AACAACCTCA	AGTACAACCTG	305	CGATACCTAG	CTGCCTTAGT
266	CGGCTGTAAC	CAAGTCTGCA	306	GCTATGACCT	GATGTGGTAC
267	GTAGCAGGTG	GTGCAGCTAC	307	TTCGCTGGTC	TCCTACAGCG
268	TCTTGTACGT	TCCTGTGAGT	308	AAGCGATAGA	AGAAGACCTA
269	AACAACCTTG	ACCATCCGAA	309	CGAGCACTGT	AGTCTAAGGA
270	CCACCACCTA	CATTTCGATTC	310	GCTATTGATC	GAGGACGATC
271	GTGGTGAGAC	GTAGGATAGG	311	TAGTACACAG	CCATGTTTCAT
272	TGTTGTGACT	TGGCATGCCT	312	ATCCGGTGCA	TTCACGCTCG
273	CACGATGTGA	AGTCATGACA	313	GTTTCATCGAG	GTACAATTTCG
274	TGGTCCAATT	TTGAGCCTGG	314	TACACATCGC	CATTCCGCAA
275	ATTAGGTGCC	CAAGTATCTC	315	ACGTGCGACT	TGGATGCAGC
276	GCACTACCAG	GCCTCGAGAT	316	CGAGTGATTA	ACCGGTAGTT
277	AAGTGATCCA	CACTCGTAGC	317	CAGCTCCTCT	ATGCGTTCCG
278	GTACATCTTC	GCTGGTACCT	318	GCAGATAAGC	CCAACACTAC
279	TCTGCCGGAT	AGGATCCTAA	319	TGTTGATGAG	GACGTCGATA
280	CGCATGAAGG	TTACAAGGTG	320	ATCACGGCTA	TGTTAGAGGT

321	CCATACTCAG	TGGTTCGTCG	353	CTTGGTGACC	ACACTAGCCT
322	GTGAGTATTC	ACCGCACGTT	354	AGAACACTTG	TGTGAGCGAA
323	TGTGCGGACA	CTTAAGTCAC	355	GCGCACACAA	CACTCCATGC
324	AACCTACGGT	GAACGTAAGA	356	TACTTGTGGT	GTGAGTTATG
325	TGAACCGAAG	CTCCGCAAGT	357	TTCAACAAGC	TACTCGGAAG
326	ACCTGGCGTA	GATGTTGCAA	358	CATGCACGCT	ATTGGCCTGC
327	CAGCTATCCT	TGATAATGCG	359	GCATGTGTTA	CGGCATTGTT
328	GTTGATATGC	ACGACGCTTC	360	AGGCTGTCAG	GCAATAACCA
329	CCGCTTGTTG	CGGACGTAA	361	TTGTACGGCT	GTAGCGAGTT
330	TTCTGCTACA	ACTCAACGGC	362	CGAGTTAAGG	AACCGTTCGA
331	AGAACAAGGT	GTATGCAATG	363	GATCCACTTA	TCTATCGTAC
332	GATGAGCCAC	TACGTTGCCT	364	ACCAGGTCAC	CGGTAACACG
333	GCTCTATCCG	GCACTTCATG	365	TCCACGGTTC	GCAACTTGTC
334	TGAGCTAGTA	CTGGACAGGT	366	ATATTCTCGG	ATGGAGACGG
335	CTGTAGGAGC	TGTTGGTCCA	367	CGGCGAAGAT	CGTCTAGTAA
336	AACAGCCTAT	AACACAGTAC	368	GATGATCACA	TACTGCCACT
337	GCGCATCGTT	TTCAAGGAGA	369	GACCTACTAT	GATCAACAAG
338	TTAAGCATCC	CCGTCTAGTG	370	TGGAATACCA	ACATTCTTCC
339	CGCGTATCAG	AGAGGATCCT	371	ATTGCCGAGC	TTGGCTGCGT
340	AATTCGGAGA	GATCTCCTAC	372	CCATGGTGTG	CGCAGGAGTA
341	AGCGAATGAT	ACCATGAGAC	373	AGATTGTGCG	CGTGTATGTC
342	CATCGTCCTG	TGTTCATCTG	374	GATGAAGTGC	GTCTGTCAAG
343	GCGTTCAAGA	CTAGACCTGT	375	TCGACTACAT	TCAACCGCCA
344	TTAACGGTCC	GAGCGTGACA	376	CTCCGCCATA	AAGCAGATGT
345	CGGATATCAA	TGCTGGATCT	377	AGTTGCCGGA	GACAAGTGGT
346	GTACCGAAGC	AAGGCACAAC	378	CCAACTTAAC	CTGGTACACA
347	TACTACCGTT	CTAATCTGGA	379	TAGGTGATCT	TGACCTGCTG
348	ACTGGTGTCG	GCTCATGCTG	380	GTCCAAGCTG	ACTTGCATAC
349	CGTAATTCAG	GATCCTCAAC	381	TCCAAGCTCC	AGTGAGGTCA
350	GAAGCGGTCA	TGCAGCAGCA	382	AGGTTGCGTG	CTACTAAGGC
351	TTCCGAAGGT	CTAGTGTTGG	383	CTAGGTAAGA	GAGTCCTATT
352	ACGTTCATC	ACGTAAGCTT	384	GATCCATCAT	TCCAGTCCAG

The dual-index sequences on MGI platform are designed to be balanced within each group. To achieve the best sequencing quality, it is recommended to use continuous indexes as indicated within one block in the table and involve at least 8 indexes for multiplexed sequencing.

